

Automated Size Selection of NEBNext® Small RNA Libraries with the Sage Pippin Prep



Daniela Munafò¹, Sadaf Hoda², Laurence Ettwiller¹, Brad Langhorst¹, Eileen Dimalanta¹, Fiona Stewart¹, Chris Boles²

Introduction

One of the fastest growing areas of biological research is regulatory small RNA structure, processing, and function. Next generation sequencing (NGS) is the method of choice for studying the variety and expression of small RNAs.

A common problem in NGS methods for small RNAs is contamination from adapter-dimer artifacts, because the artifacts are very close in size to the small RNA library elements. To address this problem, the NEBNext Small RNA Library kits from New England Biolabs use specially engineered RNA ligases and optimized workflows, and novel technology (patent pending) to dramatically reduce the formation of adapter-dimer artifacts during library construction. Since the workflow uses total RNA as the starting material, it is beneficial to perform a final size selection step on the amplified libraries. NEB has previously validated size selection methods using AMPure XP beads and manual preparative gel electrophoresis on a 6% polyacrylamide gel. Here we validate the use of Pippin Prep 3% gel cassettes for size selection of NEBNext Small RNA libraries.

Methods:

NEBNext Small RNA libraries for Illumina were constructed according to the [NEBNext Small RNA Instruction Manual](#). Following PCR amplification, the libraries were QC'd on an Agilent Bioanalyzer[®]. 1 µl of the purified PCR reaction on an Agilent Bioanalyzer using a DNA 1000 chip, according to the manufacturer's instructions. The miRNA library should appear as a peak at 147-149 bp (for a 21 nucleotide insert; Figure 1A). Size selection was performed either manually, using polyacrylamide gel electrophoresis (PAGE), or using an automated agarose gel electrophoresis on a Pippin Prep (Sage Science). The manual PAGE was performed according to the NEBNext Small RNA Instruction Manual.

Size selection using the Pippin Prep was carried out as outlined below:

Concentrate the PCR-amplified NEBNext small RNA library (100 µl)

Use a QIAQuick[®] PCR Purification Kit, with the following protocol modifications:

- Before eluting the DNA from the column, centrifuge the column with the lid of the spin column open for 5 minutes at 13,200 rpm, in order to remove all traces of ethanol from the column.
- Elute amplified DNA in 30 µl nuclease-free water.

Size selection of concentrated small RNA library using a Pippin Prep 3% DF cassette

General instructions for setting up runs on the Pippin Prep can be found in the [Operations Manual](#) (available from <http://www.sagescience.com/support>)

Create a protocol

1. In the Pippin Prep software, go to the Protocol Editor Tab.
2. Click "Cassette" folder, and select "**3% DF Marker F**". (Note: This is an internal standard cassette. Markers are provided premixed with the sample loading solution.)
3. Select the collection mode as "Range" and enter the size selection parameters as follows: BP start (110) and the BP end (160). BP Range Flag should indicate "broad".
4. Click the "Use of Internal Standards" button.
5. Make sure the "Ref Lane" values match the lane numbers.
6. Press "Save As" and name and save the protocol.

Prepare a 3% DF Pippin Prep Cassette

Follow the instructions in the Pippin Prep Operations Manual (Chapter 6) for preparing and testing a cassette.

¹ New England Biolabs, Ipswich, MA; ² Sage Science, Beverly, MA.




Figure 1: Library size distribution before and after size selection. Bioanalyzer traces from NEBNext Human Brain miRNA libraries before size selection (A) and after size selection on a 6% polyacrylamide gel (B) or a 3% agarose (dye-free) cassette for pippin prep (C). Instrument Program Mode = Range; Start (bp) = 110 and End (bp)= 160 (C).

Load and run samples on Pippin Prep

1. Bring loading solution to room temperature.
2. For each sample, combine 30 μ l of library sample with 10 μ l of DNA marker H.
3. Mix samples thoroughly (vortex mixer). Briefly centrifuge to collect.
4. Load 40 μ l (library sample plus marker) in each sample well of the 3% agarose cassette.
5. Run the program with the settings indicated above.
6. After sample has been eluted, collect size-selected sample (40 μ l) from elution module.

QC size-selected products

Run 1 μ l of the size-selected library in an Agilent Bioanalyzer using the High Sensitivity chip. (Purification of the size-selected material is not required for QC in the Bioanalyzer.)

Data Analysis

Trimmed and length-filtered reads were first mapped to miRBase (Release 20, human) using bowtie2 (option--sensitive-local --non-deterministic) and the remaining reads that do not map to miRBase were further mapped to the human genome (hg19) using Star (option--outFilterMatchNmin 20 –

outFilterMismatchNoverLmax 0.05 --outFilterMultimapNmax100000). The position of the mapped reads were compared to genomic features described in gencode annotation (v17).

Results:

Library Profile before and after size selection

The experiments were designed to compare NEBNext Small RNA library quality using two gel-based protocols: manual preparative electrophoresis on 6% PAGE gels, and automated agarose gel electrophoresis on the Pippin Prep.

Six human brain small RNA libraries were prepared from total RNA and indexed using different barcodes. Three libraries were PAGE-size-selected (**Figure 1B**) and the other three libraries were size-selected on the Pippin Prep using a 3% dye-free agarose cassette (**Figure 1C**). Both size selected libraries have only a single peak at ~149 bp, with minimal contamination from smaller or larger species. Library yields from the size-selected libraries using the manual gel and Pippin Prep was similar (~20% of the input material is recovered; **Figure 2**).




Figure 2: Library yield after size selection.

The percentage of library yield recovered after size selection is comparable for both size selection methods.

MiSeq® Sequencing and Data Analysis

The 6 barcoded-libraries were pooled in equimolar concentrations, loaded onto a MiSeq reagent v2 kit at 8 picomolar final concentration, and sequenced on MiSeq instrument (SE; 1X 36bp; 2.5 million of reads/library).

Adapter trimming and filtering reads by length

Reads were adapter-trimmed and filtered by length. Reads shorter than 15 nucleotides were discarded. A high percentage of reads (> 90% total reads for both size selection methods) passed the length filtering. The NEBNext Small RNA library prep improved workflow prevents adapter-dimer formation, therefore only a minimal percentage of reads did not contain insert.

Mapping Rate

A high percentage of trimmed reads (>66%) mapped to human miRNAs present in miRBase for both gel and Pippin Prep size selection methods (**Figure 3**). The vast majority of trimmed reads (98%) mapped either to the human genome or to miRBase. From the reads that map to human, an average of 28% overlapped at least one exonic feature (miRNA, lincRNA, pseudogene, snRNA, protein coding RNA, snoRNA, rRNA, sense_intronic, mt-tRNA, antisense RNAs) (**Figure 4**). Some of the reads that did not map to miRBase overlapped with a putative microRNA, increasing the number of reads mapping to various microRNA (from miRBase or gencode annotations) to more than 68% of total trimmed reads.




Figure 3: Mapping Rate. Libraries have high percentage of reads mapping to miRBase.




Figure 4: After aligning to miRBase, remaining reads were mapped to the human genome. The majority of the reads (72%) mapped to un-annotated regions of the genome. The remaining 28% mapped to small nucleolar RNA (snoRNA), miRNA, misc_RNA, ribosomal RNA (rRNA), small nuclear RNAs (snRNA), protein_coding RNA, Mitochondrial tRNA (Mt_tRNA) and Mitochondrial ribosomal RNA (Mt_rRNA). Very few reads mapped to pseudogenes, long intronic RNA (lincRNA), processed_transcripts and antisense RNA.

miRNA Expression Analysis

miRNA expression correlation of the 500 most abundant miRNAs between PAGE and Pippin prep size-selected libraries was excellent [$R^2 = 0.99805$] ($y = 1.0558x + 15.089$) (Figure 5). This data indicates that small RNA expression levels were not biased due to different size selection methods.




Figure 5: miRNA correlation. miRNA expression correlation was excellent between both size selection methods (N= 500 most abundant miRNAs).

Discussion and Conclusions:

The experiments reported here show that NEBNext Small RNA libraries produced using the Sage Pippin Prep equal or exceed the quality of libraries size selected using the manual gel procedure. The benefits are seen in library purity.

Two other key features of the Pippin Prep that are not shown in the present experiments are reproducibility and ease-of-use. The manual gel procedure has many individual manual steps that are time-consuming and extremely difficult to perform reproducibly. In contrast, the Pippin Prep procedure requires no manual manipulations except for gel loading. The entire separation and size-selection process is controlled by the onboard computer. This removes all opportunity for the user to introduce variability into the process.

Hands-on effort and time for the Pippin Prep procedure is 15-20 minutes to set up and load a cassette, and less than 5 minutes to remove samples at the end of the run. Run time for the NEBNext Small RNA libraries in the 3% internal standard cassette is a little over one hour. Up to five samples can be run per cassette.

In summary, we have developed an optimized Pippin Prep protocol for use with the NEBNext Small RNA library kits. The Pippin Prep protocol provides all of the sequence quality benefits of the standard manual gel protocol, but with greatly enhanced ease-of-use and reproducibility.

For more information about NEBNext Small RNA Library kits see <http://www.nebnext.com> and/or inquire by email to nebnext@neb.com

For more information about Pippin Prep see <http://www.sagescience.com> and/or inquire by email to info@sagescience.com

Sage Science, Inc.
Suite 3150
500 Cummings Center
Beverly, MA 01915-6147
info@sagescience.com
www.sagescience.com
978.922.1832

New England Biolabs, Inc.
240 County Road
Ipswich, MA 01938-2723
info@neb.com
www.neb.com
Tel: 978.927.5054
Orders: 800.632.5227
Tech: 800.632.7799
Fax: 978.921.1350

