

Guidelines for Preparing HLS-CATCH Enzyme Mix

A. Reagents: Supplied by User

The HLS-CATCH system has been tested and validated using the following reagents and suppliers:

S.pyogenes Cas9 enzyme, wild Type:

New England Biolabs		
M0386T	400 pmol	20 μM
M0386M	2,000 pmol	20 μM

Guide RNAs:

Integrated DNA Technologies (IDT)

These are provided in two halves; crRNA and tracrRNA, and must be annealed before they will assemble with the S.pyogenes Cas9 enzyme. The crRNAs must be designed to flank the target sequence, and the tracrRNA is identical for all gRNAs.

B. Preparation: 30 minutes

i The HLS-CATCH enzyme mix can be prepared prior to starting the HLS Workflow, or during the Extraction Stage of the Workflow which typically requires 1 hour. If the enzyme mix is prepared ahead of time and placed on ice, it will remain stable for several hours.

1. Anneal the tracrRNA and crRNAs

- a. Add equimolar concentrations of tracrRNA and crRNAs. Stocks are dissolved in IDT Duplex Buffer (30mM HEPES, pH7.5, 100 mM Potassium Acetate):

	vol. μl	stock [] μM	final [] in anneal mixture μM
tracrRNA	20	100	50
crRNA1	10	100	25
crRNA2	10	100	25
TOTAL	40		

- b. Heat the mix at 95oC for 5 minutes (in a heating block or thermal cycler)
- c. Remove the mix from heat, and allow to cool on the bench-top



In some cases it might be useful to use gRNAs which cut at closely spaced sites on either side of the target sequence. In this instance, anneal with equimolar amounts of each crRNA, with the total moles of crRNA equal to the total moles of tracrRNA. If designing a system with two cut sites on either side of the target, anneal with all four crRNAs at 12.5 μ M and tracrRNA at 50 μ M.

2. Assemble the Cas9-gRNA complete reaction mixture

a. Using the following order of addition, assemble the Cas9-gRNA reaction mixture:

order of addition	vol. μ l	stock [] μ M	final []	
			in enzyme mixture μ M	
1	16			Nuclease free H ₂ O
2	20	4	1	4X Enzyme Buffer (provided by Sage Science*)
3	40	50	25	Annealed gRNA mix
4	4	20	1	NEB Cas9 nuclease, wt
TOTAL	80			

b. Mix by pipetting up and down.

c. Incubate at 37°C for 10 minutes in a thermal cycler or heat block.

d. The Enzyme Mix can be placed on ice for up to 3 hours prior to use.

* provided with the HIT-0004 and HIT-0012 Cassette Kits.