

Fragments like you like them.

Pippin LINK™ Group Control



Size Select up to 105 samples per day with 3 Pippins

Pippin Link Group Control is the best way to reproducibly provide the highest quality size-selected DNA in high-throughput. By allowing operation of 1-8 Pippin Prep or BluePippin instruments from a laptop PC, bench space is conserved, and workflow is streamlined.

Sage Science offers attractively priced bundles that include multiple instruments along with the Pippin Link Control Kit. Product bundles include a laptop computer, software, and connectivity for the control of multiple units, and the Pippin instruments. The system is pre-configured and is easily user-installed. Contact us for more information or a price quote.

Ordering Information:

<u>Product Number</u>	<u>Description</u>
LNK0001	Pippin Link controller. Includes laptop*, software, and ethernet hub
PIPLINK	Pippin Prep instrument, configured for Pippin Link.
BLULINK	BluePippin instrument, configured for Pippin Link.

*also includes an acrylic stand to place the laptop over a Pippin