

PippinHT Software Upgrade Release notes

Current software version: 1.12 cassette definition set 14

Release notes:

Date: 06/01/18

Version: 1.12

v1.12-HT-CD14-060118

What's new from v.1.08 to 1.12?

- Timeout values adjusted for the following cassette definitions:
 - 3% Agarose 90-250bp 30A
 - 3% Agarose 100-250bp 30G
 - 0.75% Agarose 6-10kb high-pass 75D
 - 0.75% Agarose 15-20kb high-pass 75E
 - 0.75% Agarose 30-40kb high-pass 75F
- Bug fix for timeout with ignored peaks
- Wording change on low voltage pause

Release notes:

Date: 07/20/16

Version: 1.08

v1.08-HT-CD9-072016

What's new from v.1.07 to 1.08?

- Peak calling fail-safes added.
 - Included 0.75% Agarose 6-10kb high-pass 75D cassette definition
 - Included 0.75% Agarose 15-20kb high-pass 75E cassette definition
 - Included 0.75% Agarose 30-40kb high pass 75F cassette definition
 - Included 1.5% Agarose 300-1500bp 15C cassette definition
 - Included 1.5% Agarose 1000bp high-pass 15C cassette definition
 - Included 2% Agarose 100-600bp 20B cassette definition
 - Included 3% Agarose 90-250bp 30A cassette definition
 - Included 3% Agarose 100-250 30G cassette definition
 - Included rinse cassette definition
 - Extended delay time to eliminate initial signal spike to 15s
 - Included front LEDtest button in Factory Setup
-